Troubleshooting of infectious tests serology and immunoassay

D. mehrdad vanaki

DCLS / EMBA

Trainer & consultant QMS

Diagnostic technique	Principle	Strengths	Weaknesses	Variants	References
Immunoassay	Formation of antigen–antibody through recognition and binding	High sensitivity High specificity High speed throughput Quick TAT (20 min or less) Different types of tags, labels Automated method	Must rely on QC assurance High risk of interferences High cost	RIA EIA (FPIA, MEIA, CLIA)	Gupta <i>et al.</i> (2015), Mixson-Hayden <i>et al.</i> (2015)
NAAT	Amplification and detection of sequences from the viral genome (DNA or RNA)	High sensitivity High specificity Multiplexed platforms Genotyping Determination of the viral load New compact and portable formats	Longer run time Requires specific primers for the targets	RT-PCR qPCR NASBA TMA	García-Arroyo <i>et al.</i> (2016), Reijlans <i>et al.</i> (2008), Renois <i>et al.</i> (2010), Afshar and Mollaie (2012), Mercier-Delarue <i>et al.</i> (2014), Wu <i>et al.</i> (2014)
NGS	Polymerization of DNA template by incorporation of labelled dNTPs, and terminate the extension	High sensitivity High specificity Identification of novel genomic sequences Genotyping Accurate detection of mutations and drug resistant mutations	High cost Needs bioinformatics skills for data analysis Delay in use for routine clinical diagnostics	Pyrosequencing Fluorescently labelled dNTP Detection of released hydrogen ion (H ⁺)	Capobianchi <i>et al.</i> (2013), Bartolini <i>et al.</i> (2015), Lowe <i>et al.</i> (2016), Liu <i>et al.</i> (2013), Van den Hoecke <i>et al.</i> (2015), Rothberg <i>et al.</i> (2011)
MS	Ionization of the sample, then separation and detection of the particles according to their mass-to-charge ratio (<i>m/z</i>)	High sensitivity Versatility Cost-effectiveness High workload	Expensive equipment Limited database library	MALDI-TOF MS ESI MS Often combined with other methods: PCR-MS	Lévêque <i>et al.</i> (2014), Mengelle <i>et al.</i> (2013), Qian <i>et al.</i> (2014)

CLIA, chemiluminescent immunoassay; dNTP, deoxyribonucleotide triphosphate; EIA, enzyme immunoassay; ESI, electrospray ionization; FPIA, fluorescence polarization immunoassay; MALDI-TOF, matrix-assisted laser desorption ionization time-of-flight; MEIA, micro-particle enzyme immunoassay; MS, mass spectrometry; NGS, next-generation sequencing; NAAT, nucleic acid amplification test; NASBA, nucleic acid sequence-based amplification; qPCR, quantitative polymerase chain reaction; RIA, radio-immunoassay; RT-PCR, real-time polymerase chain reaction; TMA, transcription-mediated amplification.

PCR infectious pannel

PCR infectious pannel detect over 100 micdroorganism by qualitative and quamtitative method

PCR infectious test use foe detect slow growth organism and microorganism impossible for culture method or new strain of microorganisms

In PCR method we can use all source of clinical samples (fresh or old)


PCR advantage from culture:

Cultivation of samples for at least 24 hours and Sometimes **it takes several weeks to become positive**. If **sick If he has taken antibiotics, the culture may be negative**. Even if clinical samples are sent to the laboratory later

Cultivation of some infectious agents is negative

PCR infectious pannel

- Culture of some infectious agents such as viruses, chlamydia, Mycoplasma and rickettsiae could not be cultured by conventional methods and need specific culture media so with the help of PCR, the above samples can be checked and the result will be ready in a few hours
- An example is tuberculosis culture as gold standard, which takes time several weeks, but with the 24-hour PCR technique, the patient's result can be determined so Treatment of the patient begins quickly
- For positive TB smears, 98% PCR test is sensitive, but for negative TB smears There is 45% sensitivity and the test specificity is above 98%
- For two groups of patients PCR tuberculosis have not order , the first group of patients who have started anti-tuberculosis treatment should not be requested because the false positive PCR test occurred
- and the second group of people with negative TB smear with weak clinical suspicion of tuberculosis because the PCR result If it is false negative, the patient will not be identified.
- PCR in extrapulmonary tuberculosis is highly recommended such as CSF and serous fluid

- Detection of rising titres of antibody between acute and convalescent stages of infection,
- or the detection of IgM in primary infection.

- **Terms used in evaluating test methodology :**
- **Analytical Sensitivity** – ability of a test to detect very small amounts of a substance
- **Clinical Sensitivity** – ability of test to give positive result if patient has the disease (no false negative results)
- **Analytical Specificity** – ability of test to detect substance without interference from cross-reacting substances
- **Clinical Specificity** – ability of test to give negative result if patient does not have disease (no false positive results)

- **Affinity** : Affinity refers to the **strength of binding between a single antigenic determinant and an individual antibody combining site**. Affinity is the equilibrium constant that describes the antigen-antibody reaction
 - It is the sum of the attractive and repulsive forces operating between the antigenic determinant and the combining site .
- **Avidity** : Avidity is a measure of the **overall strength of binding of an antigen with many antigenic determinants and multivalent antibodies**
Avidity is influenced by both the valence of the antibody and the valence of the antigen.
- **Avidity is more than the sum of the individual affinities.**

Preanalytical error infectious tests

- Icteric & lipemic sample
- Clot or fibrin in samples
- Missing sample
- Mislabeling sample
- Insufficient samples
- No drug using history for patients (especially in TORCH panel)
- Dislocation of serum in sampling or separation phase

Analytical error infectious tests

- Fluctuation of temperature in refrigerator or freezer
- Frequent freeze and defreeze of serum samples
- Using expired kits and reagents
- Don't use external control
- Fluctuation of environmental factors such as moisture and temperature for manual immunoassay and ELISA methods
- Do not reach the ambient temperature of the kits taken out of the refrigerator at the time of use
- Dilution error in ELISA kits
- Error temperature in ELISA methods (37 / RT)
- Error temperature of incubator of automatic instrument for immunoassay tests (CLIA / ECL / ELFA)

Analytical error infectious tests

- Contamination and drop of conjugate reagent in manual ELISA methods
- Contamination and drop trigger reagent in automated immunoassay methods
- Carry over or well to well contamination between neg serum with pos control or between pos serum with neg serum or between neg control with pos control in ELISA methods
- Carry over in automated immunoassay methods because of fail in washing process
- Decrease or increase the number of washes in ELISA manual methods and use old and unsuitable buffer for washing
- Use of kits with improper sensitivity to test the torch and increase the amount of false negatives
- Use of kits with improper specificity to test the infectious tests and increase the amount of false positive.

Analytical error infectious tests

- Improper soaking time in manual ELISA methods
- Error in sampling and dispensing serum in manual ELISA methods
- Error in sampling and dispensing serum in automated immunoassay instruments because of complete or partial block of Probe sample or reagents with fibrin clot or Foreign object
- Lack of attention to the exact cut-off formula in ELISA manual methods, especially torch and hepatitis qualitative tests
- Displacement of serum sampling and dispensing in manual ELISA methods

Post analytical error infectious tests

- Transcriptional errors in manual data entry
- Prolonged TAT of emergency infectious tests for patients from predetermined limits by lab
- Missing or Prolonged reporting and notification of critical results in field of infectious test (especially pos viral markers for dialysis or presurgery)
- Displacement of final report during releasing results(release to wrong patients)
- Attach another mistake report to main report of patients
- Incomplete archiving or failure to maintain positive infectious specimens for repeat or legal issues
- Wrong reference interval of infectious tests in order of change of kits
- Release of pos final infectious results without reconfirm with same method or confirmatory method (if available)
- Release of false negative infectious results without check of history or deltacheck

Post analytical error infectious tests

- Lack of attention to the history of taking certain drugs for certain infectious tests such as torch (acetaminophene/ aspirin / metformin / ibuprofen)

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
Trainer & consultant QMS




MAIN STRUCTURE OF CORONAVIRUSES

Structural Protein	Function of Protein
Nucleocapsid Protein (N)	Bound to RNA genome to make up nucleocapsid
Spike Protein (S)	Critical for binding of host cell receptors to facilitate entry of host cell
Envelope Protein (E)	Interacts with M to form viral envelope
Membrane Protein (M)	- Central organiser of CoV assembly - Determines shape of viral envelope

- It has been noted that some CoVs do not need to have the full ensemble of structural proteins to make virions, highlighting that certain proteins may be dispensable or compensated by the function of non-structural proteins.

DIFFERENT TYPE OF ANALYTE LABORATORY TESTING RELATED COVID-19

RNA

ORF1a/b Gene
N Gene
E Gene
S Gene

Antigen

N Protein
S Protein

Antibody

IgM
IgG
IgA

Host Respons

CBC
CRP
D Dimer
SGOT
Albumin
LDH, etc

Potential Succability

ACE2 Gene
HLA Gene

rRTPCR, LAMP,
NGS

ELISA, Immuno-
chromatography

ELISA, Immuno-
chromatography,
Chemiluminescence
immunoassay

Enzymatic,
colorimetry,
flowcytometry,
impedance

Genotyping
microarray, RTPCR,
Sanger Seq, NGS

FAVORABLE TECHNOLOGY DETECTION RELATED SARS COV-2

Molecular Based testing

- **Closed System**
 - More Safety & Standardize
- **Open system**
 - Easy & Faster to develop new test

Immunoassay Based testing

- **Rapid immunochromatigraphy**
 - Faster results, Easy to use
- **ELISA or Chemiluminescence**
 - More standardize & possible to Quantify

TESTS FOR SARS–COV-2/COVID-19 AND POTENTIAL USES

Report from the American Society for Microbiology COVID-19 International Summit, 23 March 2020: Value of Diagnostic Testing for SARS–CoV-2/COVID-19

Type of Test	Measure	Value	Beneficiary
 Nucleic acid amplification test for viral RNA <i>(nasopharyngeal swab, oropharyngeal swab, sputum, bronchoalveolar lavage fluid, others)</i>	Current infection with SARS-CoV-2	• Inform individual of infection status so they can anticipate course of illness and take action to prevent transmission • Inform patient management and actions needed to prevent transmission • Inform actions needed to prevent transmission	• Individual • Healthcare or long-term care facility • Public health
 Antibody detection	Past exposure to SARS-CoV-2	• Detect susceptible individuals (antibody negative) and those previously infected • Identify individuals with neutralizing antibodies • Facilitate contact tracing and surveillance	• Identify those potentially immune to SARS-CoV-2 (if tests can detect protective immunity, individuals could be returned to work) • Healthcare facilities: Experimental therapy • Public health

RECOMMENDED TEST LIST:

Laboratory Test	Main laboratory abnormalities observed in adult patients with unfavorable COVID-19 progression (Modified from 1-12)	Potential clinical and biological significance (Modified from 3)
Complete blood count	Increased white blood cell	Bacterial (super)infection
	Increase neutrophil count	Bacterial (super)infection
	Decreased lymphocyte count	Decreased immunological response to the virus
	Decreased platelet count	Consumption (disseminated) coagulopathy
Albumin	Decreased	Impairment of liver function
Lactate Dehydrogenase	Increased	Pulmonary injury and/or widespread organ damage
Alanine Aminotransferase	Increased	Liver injury and/or widespread organ damage
Aspartate aminotransferase	Increased	Liver injury and/or widespread organ damage
Total bilirubin	Increased	Liver injury
Creatinine	Increased	Kidney injury
Cardiac troponin	Increased	Cardiac injury
D-Dimer	Increased	Activation of blood coagulation and/or disseminated coagulopathy
Prothrombin Time	Increased	Activation of blood coagulation and/or disseminated coagulopathy
Procalcitonin	Increased	Bacterial (super)infection
C-reactive protein	Increased	Severe viral infection/viremia/viral sepsis
Ferritin	Increased	Severe inflammation
Cytokines (IL-6)	Increased	Cytokine storm syndrome

FREQUENT LABORATORY ABNORMALITIES IN PATIENTS WITH COVID-19 INCLUDE:

- **Lymphopenia** (35-75%)
- **↑ C reactive protein** (CRP **75-93%**)
- **↑ Lactate dehydrogenase** (LDH; 27-92%)
- **↑ Erythrocyte sedimentation rate** (ESR; up to 85%)
- **↑ D-dimer** (36-43%)


- **↓ Albumin** (50-98%)
- **↓ Hemoglobin** (41-50%)

SEROLOGICAL TESTING FOR SARS CORONA 2

- There has been **much debate regarding the current value** of serological testing in COVID-19 **diagnosis and monitoring**.
- Serologic based tests are **not currently recommended** by the **CDC, NHS or other health organizations**.
- There is general concern regarding their use in the acute phase of infection as they detect infection too late in the course of illness (usually more than 7-10 days) and they also may cross-react with serologic responses to seasonal coronaviruses.
- However, there is anticipated value in using improved serological testing in the future **for public and occupational health monitoring and assessment**

Timing of Antibody Response to SARS-CoV-2

- New virus = no pre-existing antibodies or immunity
- We are still learning about our immune response to SARS-CoV-2
 - Many develop Abs ~1-2 weeks after symptoms
 - Due to delay in seroconversion, Abs do not play a routine role in diagnosis
 - >95% of patients are Ab positive after 2 weeks
 - Some patients may not seroconvert
 - Immunostatus
 - Assay dependent?
 - Severity of illness?
 - IgM declines 5-7 weeks post onset
 - IgG remains positive for ≥ 10 weeks post onset

SEROLOGICAL TESTING FOR SARS CORONA 2

- For **serum antibody monitoring**, according to the “**Handbook of Covid-19 Prevention and Treatment**” (Zhejiang University School of Medicine)
- **Serum IgM** is detectable **10 days** after symptom onset
- **Serum IgG** is detectable **12 days** after symptom onset.
- A positive interpretation has been defined as a positive IgM, or an increased IgG titer > 4 times than that in the acute phase.

Rt-PCR is Considered the “Gold Standard”

- **98% sensitivity and 99% selectivity**
- **Requires 20 minutes to 90 minutes to run the test**
- **Turnaround: 24 hour to >7 days; very problematic**
- **Three possible results:**
 - **detected;**
 - **inconclusive**
 - **not detected**

NUCLEIC ACID TESTING (RNA SARS CORONA 2)

Upper airway specimens

Nasopharyngeal swabs

Nasal swabs

Nasopharyngeal secretions

Lower airway specimens

Sputum

Airway secretion

Bronchoalveolar lavage fluid

Other

Blood

Feces

Urine


Conjunctival secretions

collected by a healthcare professional, If both NP and OP swabs both are collected, they should be combined in a single tube to maximize test sensitivity and limit testing resources

Gold standard for Diagnosis COVID-19

Real-time reverse transcription polymerase chain reaction (rRT-PCR) is the current gold standard for diagnosing suspected cases of COVID-19

Institute	Gene targets
China CDC, China	ORF1ab and N
Institut Pasteur, Paris, France	Two targets in RdRP
US CDC, USA	Two targets in N gene (previously 3)
National Institute of Infectious Diseases, Japan	Pancorona and multiple targets, Spike protein
Charité, Germany	RdRP, E, N
HKU, Hong Kong SAR	ORF1b-nsp14, N
National Institute of Health, Thailand	N

COMPARISON WHO VS. CDC METHOD FOR RTPCR SARC CORONA 2

Test	Molecular targets	Scope	Limit of blank	Reference specimens	Storage conditions
WHO					
	E gene	First-line screening	3.9 copies × reaction	Nasopharyngeal AND oropharyngeal swab or wash in ambulatory patients, lower respiratory specimens (sputum and/or endotracheal aspirate or bronchoalveolar lavage)	≤5 days: 2–8 °C
	RdRp gene	Confirmatory testing	3.6 copies × reaction		>5 days: ≤70 °C
	N gene	Additional confirmatory testing	N/A		(dry ice)
CDC					
	N1/2/3 gene	Combined assay	1.0–3.2 copies/μL	Nasopharyngeal AND oropharyngeal swabs, sputum, lower respiratory tract aspirates, bronchoalveolar lavage and nasopharyngeal wash/aspirate or nasal aspirate	≤4 days: 4 °C
	RNase P gene	Control assay	N/A		>4 days: ≤70 °C

E gene, envelop gene; N gene, nucleocapside gene; RdRp gene, RNA-dependent RNA polymerase gene; RNase P gene, human RNase P gene.

Troubleshooting and variables

**covid -19
lab tests**

**D. mehrdad vanaki
DCLS / EMBA
Trainer & consultant QMS**

WHY MIGHT COVID-19 TESTS FAIL?

They may be in the **early stage** of the disease with a **viral load that is too low to be detected**.

They may have **no major respiratory symptoms**, so there could be **little detectable virus** in the patient's throat and nose.

There may have been a **problem with sample collection**, meaning there was very little sample to test.

There may have been **poor handling and shipping** of samples and test materials.

There may have been **technical issues inherent in the test**, e.g. virus mutation.

preanalytical variables

in

covid 19

Lab tests

PREANALYTICAL

General

- Lack of identification/misidentification
- Inadequate procedures for specimen (e.g. swab) collection, handling, transport

Storage

- Collection of inappropriate or inadequate material for quality or volume
- Presence of interfering substances
- Manual (pipetting) errors

Specific

- Sample contamination
- Testing in patients receiving antiretroviral therapy

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
Specimen Collection and Transport Devices in molecular diagnostic tests

- in general, the anticoagulants EDTA and acid citrate dextrose (ACD) should be used because they are most compatible with molecular methods.
- **Heparin tubes should not be used for amplification methods because heparin is a potent inhibitor of polymerase enzymes.**
- **Other specimens, such as tissues, urine, body fluids, and stool, may be transported in clean containers to the laboratory.**
- **Specimens in fixatives such as alcohols or formalin should be extracted as quickly as possible after receipt by the laboratory. It can be more difficult to obtain good-quality DNA from fixed tissues**

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Extraction and Purification of Nucleic Acids

- Nucleic acid extraction and purification procedures can be a source of errors in diagnostic testing.
- If intact nucleic acids are not isolated, false-negative results can occur.
- In addition, the extraction and purification procedures can be a source of interfering substances if reagents such as phenol or alcohols are carried over into the final eluted material containing the target nucleic acid.

Extraction and Purification of Nucleic Acids

- Nucleic acid preparation methods are a critical component of molecular testing because for most assays, high-quality DNA or RNA is critical to assay performance; maintained (especially important for RNA assays), and inhibitory substances must be eliminated
- The quality of purified DNA and RNA, as well as the concentration and purity, affect results of hybridization, amplification, and detection methods. The target nucleic acid must be intact, in a sufficient concentration to be detected by the assay, and be free of inhibitory substances.

Inhibition in Amplification Assays

- Inhibition of amplification reactions is a well recognized phenomenon in the diagnostic.
- Many substances from a variety of sources can inhibit PCR reactions. Inhibitors can be present in the clinical specimen (e.g., hemoglobin and.....) **Hemoglobin is a potent inhibitor of PCR amplification and a well-characterized inhibitor of the Taq polymerase**
- **Inhibitors are especially problematic in quantitative assays because they may cause aberrant results (inaccurate quantitation) or they may inhibit amplification altogether**

Mechanisms by which Inhibitors Function

- There are several mechanisms by which inhibitors function, usually targeting the enzymes involved in amplification or interfering with binding of the enzyme to the target nucleic acid:
- **1. Blocking active site of enzymes:** Some substances block active sites of enzymes used in amplification reactions. Heme compounds from blood function to inhibit Taq polymerase in this manner.
- **2. Inactivating enzymes:** Inhibitors can also act by denaturing or degrading the enzymes used in amplification (e.g., if proteinases are present in the sample).
- **3. Binding target nucleic acids,** preventing the target from being accessible to primers or enzymes.
- **4. Binding divalent cations:** Inhibition can occur when divalent cations (Mg^{2+}) necessary as co-factors for polymerase enzymes are bound by EDTA.

TABLE 21.1 Common PCR Inhibitors

General Source	Inhibitory Substance
Specimens	Blood—hemoglobin, IgG, lactoferrin Blood—DNA binding proteins Muscles—myoglobin Sputum ^a Stool ^a —bile salts, complex polysaccharides Tissues—melanin Urine ^a —urea
Specimen collection tubes	EDTA Heparin
Nucleic acid extraction	Detergents (Sarkosyl, SDS) Proteinases Protein disrupting agents Guanidinium salts Phenol Alcohols

^aFor many complex clinical specimens, inhibitors have not been well-characterized.

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covid 19

Lab tests

POTENTIAL PREANALYTICAL AND ANALYTICAL VULNERABILITIES IN THE LABORATORY DIAGNOSIS OF CORONAVIRUS DISEASE 2019 (COVID-19) USING (REAL TIME) RT-PCR.

Testing carried out
outside of the
diagnostic window

Active viral
recombination

Use of non-
adequately
validated assays

Lack of
harmonization of
primers and
probes

Instrument
malfunctioning

Insufficient or
inadequate
material

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Verification Requirements

CAP treats EUA assays similar to FDA cleared assays

Test Method Verification (COM.40300/COM.40325)

- **Analytical Interferences**
- **Precision**
- **Reportable Range**
- **Accuracy**

Precision

- **Closeness of agreement between independent test measures**
 - “Reproducibility/repeatability”
- **Typical sources of imprecision include differences in timing, temperature, mixing, pipetting, etc.**

Two aspects should be tested

- **Intra-assay precision**
 - Measurements collected under very similar conditions (i.e. same run)
- **Inter-assay precision**
 - Measurements collected under very different conditions (i.e. different operators, different instruments, different days, etc.)

Ideal to test concentrations at or near the level of detection

Accuracy

- **Extent to which a particular test is in agreement with a reference method or comparator**
 - “Trueness”
- **Ideal “comparator”: Specimens from patients with known COVID-19 infection (established through molecular testing)**
 - With the increasing prevalence of COVID-19 infections, most laboratories should be able to obtain these
- **Secondary “comparator”: Specimens with known positive and negative antibody status tested using another validated/verified antibody test**

Determination of Accuracy

- Qualitative test data typically analyzed in a 2x2 table

Specimen type	Reference Method		
		Positive	Negative
Method being evaluated	Positive	True positives	False positives
	Negative	False negatives	True negatives

Sensitivity: $TP / (TP + FN) \times 100$

Specificity: $TN / (TN + FP) \times 100$

Set goals for each prior to experiments (typically 95%)

Laboratory Determines:

-Reference method

-How many positive and negative samples to include

-Thresholds for acceptable sensitivity and specificity

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